

Anti-GRIN2B antibody (C-Term) (STJA0003716)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 IHC 1:1000 IP 3 μ l per 200 μ g lysate
Formulation	Lyophilized
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	24410
Gene Symbol	Grin2b
UniProt ID	NMDE2_RAT
Immunogen Region	C-Term

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining iPSC-derived neurons with Anti-NMDA Receptor NR2B green and TuJ1 red STJA0003716



Western blot NR2B subunit 180 kDa immunolabeling in rat hippocampal lysate STJA0003716