

Anti-OLIG1 antibody (STJA0003731)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was protein a purified.
Dilution Range	WB 1:3000 IHC 1:100 ICC 1:100 IP 1:50
Formulation	100 μ l in 100 mM Glycine, 200 mM Tris, pH 8.0
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	50914
Gene Symbol	Olig1
UniProt ID	OLIG1_MOUSE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot 27 kDa Olig1 protein immunolabeling in neonatal rat brain lysate STJA0003731



IHC staining rat oligodendrocytes with anti-Olig1 antibody shown in red STJA0003731