

Anti-Phospho-MAPK14-Thr180/Tyr182 antibody (STJA0003734)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Rabbit
Reactivity	Human/Mouse/Bovine/Canine/Chicken/Non-Human Primates/Rat/Zebrafish

PRODUCT PROPERTIES



Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 IHC 1:250
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	Mapk14 Uniprot ID= MK14_RAT "
-------------	---------------------------------------------------------------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot phosphorylated p38 MAPK at Thr180-Tyr182 in UV-treated HeLa cells antibody phosphorylated p38 red STJA0003734
 IHC staining human breast cancer tissue with antibody phosphorylated p38 red STJA0003734