

Anti-Phospho-PGR-Ser294 antibody [608] (STJA0003756)

GENERAL INFORMATION

| | |
|---------------|--------------------------|
| Product Type | Primary antibodies |
| Applications | WB/IHC |
| Host / Source | Mouse |
| Reactivity | Human/Non-Human Primates |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | This antibody was protein g purified culture from supernatant. |
| Dilution Range | WB 1:1000 IHC 1:1000 |
| Formulation | 100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol. |
| Isotype | IgG1 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 5241 |
| Gene Symbol | PGR |
| UniProt ID | PRGR_HUMAN |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot PR-A and PR-B isoforms phosphorylated progesterone receptor in T47D cells STJA0003756