

Anti-RXRG antibody (Hinge region) [1373] (STJA0003769)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB
	Host / Source	Mouse
	Reactivity	Human/Rat/Bovine/Mouse/Non-Human Primates
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	This antibody was protein g purified culture from supernatant.
	Dilution Range	WB 1:1000
	Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
	Isotype	IgG1
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	6258
	Gene Symbol	RXRG
	UniProt ID	RXRG_HUMAN
	Immunogen Region	Hinge region
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.