

Anti-RHO antibody [1D4] (STJA0003770)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Mouse
Reactivity	Amphibians/Most Mammals

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified culture from supernatant.
Dilution Range	WB 1:4000 IHC 1:100-1:200
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining rhodopsin protein immunolabeling in mouse retinal rod spherules STJA0003770



Western blot 39 kDa rhodopsin protein labeling in rat retina STJA0003770