

Anti-SF3B4 antibody (Full-Length) [3A1] (STJA0003775)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified culture from supernatant.
Dilution Range	WB 1:2000 IHC 1:1000 ICC 1:1000
Formulation	100 ul in PBS + 50% Glycerol and 5 mm Sodium Azide
Isotype	IgG2b
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10262
Gene Symbol	SF3B4
UniProt ID	SF3B4_HUMAN
Immunogen Region	Full-Length

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining HeLa cells using anti-SAP49 antibody red fluorescence with DAPI counterstain STJA0003775



ICC staining HeLa cells using antibody nuclear SAP49 red and vimentin green STJA0003775



Western blot 49 kDa SAP49 protein immunolabeling in HeLa cell lysate STJA0003775