

Anti-SYN1 antibody (STJA0003781)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/ICC/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 ICC 1:1000-1:20001 μ g per 200 μ g lysate
Formulation	Lyophilized
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot 78 kDa synapsin I doublet protein in rat hippocampal lysate STJA0003781



ICC staining mouse caudate neurons with synapsin antibody punctate green signal STJA0003781



ICC staining rat cortical neurons with synapsin antibody punctate red distribution STJA0003781