

Anti-Phospho-SYT1-Thr202 antibody (STJA0003791)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Rat/Bovine/Canine/Chicken/Human/Non-Human Primates/Mouse/Zebrafish

PRODUCT PROPERTIES



Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 ICC 1:400
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25716
Gene Symbol	Syt1
UniProt ID	SYT1_RAT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining rat cortical neurons phosphorylated synaptotagmin using antibody STJA0003791

Western blot in synaptotagmin rat cortical lysate phosphorylation STJA0003791