

## Anti-Phospho-SYT1-Ser309 antibody (STJA0003792)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Rat/Bovine/Canine/Chicken/Human/Non-Human Primates/Mouse/Zebrafish

### PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 ICC 1:400
Formulation	100 $\mu$ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 $\mu$ g per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">25716</a>
Gene Symbol	<a href="#">Syt1</a>
UniProt ID	<a href="#">SYT1_RAT</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining rat cortical neurons phospho-synaptotagmin Ser309 and PSD95 STJA0003792



Western blot synaptotagmin phospho-Thr202 immunolabeling in rat cortical lysate STJA0003792