

Anti-TH antibody (STJA0003817)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Sheep
Reactivity	All Mammalian/Chicken

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 IHC 1:1000 ICC 1:1000
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION



Gene ID	25085
Gene Symbol	Th
UniProt ID	TY3H_RAT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining mouse locus coeruleus tyrosine hydroxylase labeled with antibody red in rat striatal lysate STJA0003817



IHC staining rabbit retina with tyrosine hydroxylase labeled in green STJA0003817