

Anti-Vimentin antibody (STJA0003823)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Chicken
Reactivity	Human/Mouse/Rat/Bovine/Canine/Feline/Goat/Guinea Pig/Hamster/Horse/Non-Human Primates/Rabbit/Sheep/Vole

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was total igy fraction.
Dilution Range	WB 1:1000 IHC 1:500-1:5000 ICC 1:500-1:10,000
Formulation	Total IgY fraction in PBS + 10 mM Sodium Azide.
Isotype	IgY
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7431
Gene Symbol	VIM
UniProt ID	VIME_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




IHC staining mouse tumor tissue using Anti-Vimentin antibody green with E-cadherin co-stain STJA0003823




Western blot 50 kDa Vimentin protein immunolabeling in NIH 3T3 cell lysate STJA0003823




ICC staining HeLa cells using anti-vimentin antibody intermediate filament network STJA0003823




ICC staining neuron-glia cultures with anti-Vimentin green and anti-GFAP red antibodies STJA0003823




ICC staining rat neurons and glia with anti-vimentin antibody filamentous green staining STJA0003823




ICC staining fibroblasts and PC12 cells with anti-vimentin green and anti-peripherin red STJA0003823




ICC staining HeLa cells using rabbit anti-vimentin antibody intermediate filaments in green STJA0003823




ICC staining HeLa cells using rabbit anti-vimentin antibody red fluorescence STJA0003823



ICC staining HeLa cells with Anti-Vimentin antibody cytoplasmic intermediate filaments STJA0003823



ICC staining HeLa cells using anti-vimentin antibody green with DAPI nuclear stain blue STJA0003823



ICC staining HeLa cells with vimentin antibody green and SAP49 red with DAPI STJA0003823