

Anti-ARG1-ARG1 antibody (STJA0003835)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Chicken
Reactivity	Guinea Pig/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was total igy fraction.
Dilution Range	WB 1:1000 IHC 1:500
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgY
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	29221
Gene Symbol	Arg1
UniProt ID	ARG1_RAT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining mouse retina arginase greenCD206 redand DAPI nuclear stain STJA0003835



Western blot arginase I protein bands at 35 and 38 kDa in mouse liver lysate STJA0003835