

Anti-ChAT-Choline Acetyltransferase antibody (STJA0003837)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Goat
Reactivity	Chicken/Guinea Pig/Human/Mouse/Non-Human Primates/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from serum.
Dilution Range	WB 1:100 IHC 1:100-1:1000
Formulation	100 μ l in PBS + 5 mg/ml BSA + 0.2% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1103
Gene Symbol	CHAT
UniProt ID	CLAT_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot 70 kDa ChAT immunolabeling in rat hippocampal lysate STJA0003837



IHC staining rat brainstem using anti-ChAT antibody specific labeling STJA0003837