

Anti-MAP2-MAP2 antibody (STJA0003848)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Goat
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified from serum.
Dilution Range	WB 1:1000 IHC 1:1000-1:2000 ICC 1:1000-1:2000
Formulation	100 ul in PBS + 50% Glycerol and 5 mm Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4133
Gene Symbol	MAP2
UniProt ID	MTAP2_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




IHC staining rat brain stem section with MAP2 antibody neuronal perikarya and dendrites STJA0003848




Western blot MAP2 protein 280 kDa immunolabeling in rat cortical lysate STJA0003848




ICC staining MAP2 in primary mouse neurons using antibody red with Hoechst nuclear stain STJA0003848



MAP2 antibody IHC staining in hippocampal CA1 neurons from 3xTg mouse brain tissue STJA0003848



IHC staining with MAP2 magenta REST green and DAPI in mouse hippocampus CA1 STJA0003848



ICC staining with MAP2 antibody red and REST green in 3xTg primary cortical neurons STJA0003848