

Anti-MARCKS antibody (STJA0003849)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Chicken
Reactivity	Human/Non-Human Primates

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was total igy fraction.
Dilution Range	WB 1:10000 IHC 1:500-1:1000 ICC 1:500-1:5000
Formulation	Total IgY fraction in PBS + 10 mM Sodium Azide.
Isotype	IgY
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4082
Gene Symbol	MARCKS
UniProt ID	MARCS_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot HeLa lysate 80 kDa MARCKS protein immunolabeling STJA0003849