

Anti-SARS-CoV-2-S-Protein ACE 2 Binding Domain antibody [2G1] (STJA0003861)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purified culture supernatant
Dilution Range	WB 1:2000 ICC 1:1000
Formulation	100 ul in PBS + 50% Glycerol and 5 mm Sodium Azide
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot GFP-tagged SARS-CoV-2 in green
no CoV-2 staining in untransfected cells
STJA0003861



Immunofluorescence SARS-CoV-2 transfected
HEK293 cells in green and nuclei in blue
STJA0003861