

Anti-SOX2 antibody (C-Term) (STJA0003864)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Mouse/Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen affinity purified from pooled serum
Dilution Range	WB 1:1000 ICC 1:20-1:200
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6657
Gene Symbol	SOX2
UniProt ID	SOX2_HUMAN
Immunogen Region	C-Term

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot SOX2 protein immunolabeling at 35 kDa in NIH 3T3 cell lysate STJA0003864



ICC staining FLAG-tagged SOX2 in COS cells using antibody with DAPI nuclear counterstain STJA0003864