

Anti-CD71 antibody [DF1513] (STJA0003869)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB/IHC |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein G purified culture supernatant. |
| Dilution Range | WB 1:1000 IHC 1:100 |
| Formulation | 100 ul in PBS |
| Isotype | IgG1 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 7037 |
| Gene Symbol | TFRC |
| UniProt ID | TFR1_HUMAN |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



antibody WB analysis 98 kDa CD71 protein band in UV treated HeLa cell lysate STJA0003869