

## Anti-Vimentin antibody [LN6] (STJA0003872)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein L purified culture supernatant
Dilution Range	WB 1:1000 IHC 1:100
Formulation	100 ul in PBS
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7431</a>
Gene Symbol	<a href="#">VIM</a>
UniProt ID	<a href="#">VIME_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis detecting 50 kDa vimentin protein in HeLa cell lysate STJA0003872