

## Anti-Ly-6G/Ly-6C (Gr-1) antibody [RB6-8C5] {ER780} (STJA0003921)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	ER780
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">100041546</a> <a href="#">17067</a> <a href="#">546644</a>
Gene Symbol	<a href="#">Ly6c2</a> <a href="#">Ly6c1</a> <a href="#">LY6C2_MOUSE</a> <a href="#">LY6C1_MOUSE</a> <a href="#">LY6G_MOUSE</a> <

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---