

**Anti-CD90/CD90.1 antibody [OX-7] {PE/Cyanine7} (STJA0003959)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Rat                |

**PRODUCT PROPERTIES**

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.2 mg/mL  |
| Conjugation         | PE/Cyanine7  |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG1k  |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

**TARGET INFORMATION**

|             |                          |
|-------------|--------------------------|
| Gene ID     | <a href="#">24832</a>    |
| Gene Symbol | <a href="#">Thy1</a>     |
| UniProt ID  | <a href="#">THY1_RAT</a> |

**ADDITIONAL INFORMATION**

|      |   |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|