

Anti-KLRG-1 antibody [2F1] {FITC} (STJA0004064)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Syrian Hamster
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	10219 50928
Gene Symbol	KLRG1 Klrg1
UniProt ID	KLRG1_HUMAN KLRG1_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---