

Anti-KLRG-1 antibody [2F1] {EV450} (STJA0004068)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Syrian Hamster
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	EV450
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	50928 10219
Gene Symbol	Klrg1 KLRG1
UniProt ID	KLRG1_MOUSE KLRG1_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---