

Anti-CD28 antibody [37.51] {AF488} (STJA0004082)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Syrian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	AF488
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12487
Gene Symbol	Cd28
UniProt ID	CD28_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---