

## Anti-CD48 antibody [HM48-1] {PerCP} (STJA0004101)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Armenian Hamster   |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.2 mg/mL  |
| Conjugation         | PerCP  |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG  |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">12506</a>      |
| Gene Symbol | <a href="#">Cd48</a>       |
| UniProt ID  | <a href="#">CD48_MOUSE</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|