

Anti-CD106 antibody [M/K-2.7] {PerCP/Cyanine5.5} (STJA0004133)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene Symbol	Vcam1 Uniprot ID= VCAM1_MOUSE "
-------------	---

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---