

Anti-CD71 antibody [R17 217.1.3/TIB-219] {PerCP/Cyanine5.5} (STJA0004135)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC
	Host / Source	Rat
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	0.2 mg/mL
	Conjugation	PerCP/Cyanine5.5
	Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
	Isotype	IgG2ak
	Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
TARGET INFORMATION		
	Gene ID	22042
	Gene Symbol	Tfrc
	UniProt ID	TFR1_MOUSE
ADDITIONAL INFORMATION		
	Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.