

Anti-CD335 antibody [29A1.4] {PE/TR} (STJA0004261)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC
	Host / Source	Rat
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	0.2 mg/mL
	Conjugation	PE/TR
	Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
	Isotype	IgG2ak
	Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
TARGET INFORMATION		
	Gene ID	17086
	Gene Symbol	Ncr1
	UniProt ID	NCTR1_MOUSE
ADDITIONAL INFORMATION		
	Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.