

Anti-KBTBD11 antibody (STJA0004342)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2.88 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:1000-1:5000
Formulation	PBS with 0.05% NaN3 and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9920
Gene Symbol	KBTBD11
UniProt ID	KBTBB_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---