

## Anti-ARIH2 antibody (STJA0004349)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Polyclonal  |
| Concentration       | 0.6 mg/mL   |
| Conjugation         | Unconjugated  |
| Purification        | Affinity purification   |
| Dilution Range      | WB 1:500-1:2000   |
| Formulation         | PBS with 0.05% sodium azide and 50% glycerol, PH7.4             |
| Isotype             | IgG   |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">10425</a>      |
| Gene Symbol | <a href="#">ARIH2</a>      |
| UniProt ID  | <a href="#">ARI2_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|