

## Anti-TSPAN13 antibody (STJA0004354)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.8 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:1000, IHC 1:25-1:100
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">27075</a>
Gene Symbol	<a href="#">TSPAN13</a>
UniProt ID	<a href="#">TSN13_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---