

Anti-ITPR2 antibody (STJA0004425)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.9 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at-20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3709
Gene Symbol	ITPR2
UniProt ID	ITPR2_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------