

## Anti-KCTD16 antibody (STJA0004501)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/ELISA       |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES

|                     |                                                                 |
|---------------------|-----------------------------------------------------------------|
| Clonality           | Polyclonal                                                      |
| Concentration       | 0.8 mg/mL                                                       |
| Conjugation         | Unconjugated                                                    |
| Purification        | Antigen affinity purification                                   |
| Dilution Range      | WB 1:500-1:2000, IHC 1:25-1:100                                 |
| Formulation         | PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4         |
| Isotype             | IgG                                                             |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">57528</a>       |
| Gene Symbol | <a href="#">KCTD16</a>      |
| UniProt ID  | <a href="#">KCD16_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|