

Anti-HHIPL2 antibody (STJA0004518)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.9 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79802
Gene Symbol	HHIPL2
UniProt ID	HIPL2_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---