

Anti-SLC35F6 antibody (STJA0004592)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.9 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:100-1:300
Formulation	PBS with 0.05% NaN3 and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54978
Gene Symbol	SLC35F6
UniProt ID	S35F6_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---