

Anti-ATAD1 antibody (STJA0004608)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB/IHC/ELISA |
| Host / Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Polyclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Antigen affinity purification |
| Dilution Range | WB 1:500-1:2000, IHC 1:50-1:300 |
| Formulation | PBS with 0.05% NaN3 and 40% Glycerol, pH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 84896 |
| Gene Symbol | ATAD1 |
| UniProt ID | ATAD1_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|