

Anti-CARD9 antibody (STJA0004761)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC/ELISA |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Concentration | 0.6 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IHC 1:50-1:200 |
| Formulation | PBS with 0.05% sodium azide and 50% glycerol, PH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at-20°C for the long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 64170 |
| Gene Symbol | CARD9 |
| UniProt ID | CARD9_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|