



Anti-RETREG1 antibody (STJA0004770)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1.02 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54463
Gene Symbol	RETREG1
UniProt ID	RETR1_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---