

## Anti-KLK9 antibody (STJA0004872)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC/ELISA          |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                 |
|---------------------|-----------------------------------------------------------------|
| Clonality           | Polyclonal                                                      |
| Concentration       | 1.5 mg/mL                                                       |
| Conjugation         | Unconjugated                                                    |
| Purification        | Antigen affinity purification                                   |
| Dilution Range      | IHC 1:40-1:200                                                  |
| Formulation         | PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4         |
| Isotype             | IgG                                                             |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">284366</a>     |
| Gene Symbol | <a href="#">KLK9</a>       |
| UniProt ID  | <a href="#">KLK9_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|