

## Anti-BOLA1 antibody (STJA0004891)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC/ELISA          |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Polyclonal  |
| Concentration       | 0.7 mg/mL   |
| Conjugation         | Unconjugated  |
| Purification        | Antigen affinity purification                                   |
| Dilution Range      | IHC 1:40-1:200  |
| Formulation         | PBS with 0.05% NaN3 and 40% Glycerol, pH7.4                     |
| Isotype             | IgG   |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">51027</a>       |
| Gene Symbol | <a href="#">BOLA1</a>       |
| UniProt ID  | <a href="#">BOLA1_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|