

Anti-CNP antibody (Full-Length) (STJA0004898)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Goat
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:5000 IHC 1:3000
Formulation	PBS + 50% glycerol and 5 mM NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1267
Gene Symbol	CNP
UniProt ID	CN37_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Full-Length
Specificity	Specific for endogenous levels of the ~46 kDa CNP protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining with CNP in mouse striatum myelinating oligodendrocytes STJA0004898



antibody WB analysis 46 kDa CNP protein in rat hippocampal lysate STJA0004898