

Anti-SARM1 antibody (Full-Length) (STJA0004904)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES



Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified from Pooled Serum
Dilution Range	WB 1:1000 ICC 1:50
Formulation	10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 µg per ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23098
Gene Symbol	SARM1
UniProt ID	SARM1_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Full-Length
Specificity	Specific for endogenous levels of the ~79 kDa SARM1 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining SARM1 in transfected COS-7 cells using antibody with DAPI nuclear counterstain
 antibody WB analysis 79 kDa SARM1 protein detection in HeLa cell lysate STJA0004904
 STJA0004904