

Anti-IL-10 antibody [SM4919] (STJA0004919)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (R) 1:1000-1:2000/1:500-1:1000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG1k
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION


Gene ID	16153
Gene Symbol	Il10
UniProt ID	IL10_MOUSE

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Immunohistochemistry analysis of paraffin-embedded rat lung using IL-10 at dilution of 1:2000.




Immunohistochemistry analysis of paraffin-embedded rat pneumonia using IL-10 at dilution of 1:1000.



Immunofluorescence analysis of paraffin-embedded mouse inflammatory lung using IL-10 at dilution of 1:300.



Immunofluorescence analysis of paraffin-embedded rat inflammatory lung using IL-10 at dilution of 1:600.



Immunofluorescence analysis of paraffin-embedded rat inflammatory spleen using IL-10 at dilution of 1:600.