

Anti-Galectin 3 antibody [SM4925] (STJA0004925)

GENERAL INFORMATION

	Product Type	Primary antibodies
	Applications	WB/IHC/IF
	Host / Source	Mouse
	Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, R) 1:500-1:1000 IHC/IF (H, M, R) 1:2000-1:4000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION



Gene Symbol [Lgals3](#) Uniprot ID=[LEG3_MOUSE](#)"

ADDITIONAL INFORMATION


Note: STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.





Western blot analysis of Galectin-3 at dilution of 1:1000



Immunohistochemistry of paraffin embedded mouse hepatitis using Galectin-3 at dilution of 1:2000



Immunohistochemistry of paraffin embedded mouse pneumonia using Galectin-3 at dilution of 1:2000



Immunohistochemistry of paraffin embedded rat pneumonia using Galectin-3 at dilution of 1:2000