

Anti-GRP78 BiP antibody [SM4944] (STJA0004944)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:500-1:1000 IHC/IF (H, M, R) 1:1000-1:2000/1:1000-1:2000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG2a
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION



Gene ID	14828
Gene Symbol	Hspa5
UniProt ID	BIP_MOUSE

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.





Western blot analysis of BiP at dilution of 1:1000





Immunohistochemistry analysis of paraffin-embedded mouse brain using BiP at dilution of 1:2000




Immunohistochemistry analysis of paraffin-embedded rat kidney using BiP at dilution of 1:2000




Immunofluorescent analysis of paraformaldehyde-fixed mouse kidney using BiP at dilution of 1:2000




Immunofluorescent analysis of paraformaldehyde-fixed rat kidney using BiP at dilution of 1:2000



Immunofluorescent analysis of paraformaldehyde-fixed mouse brain using BiP at dilution of 1:2000



Immunofluorescent analysis of paraformaldehyde-fixed rat brain using BiP at dilution of 1:2000



Immunofluorescent analysis of paraformaldehyde-fixed rat colon using BiP at dilution of 1:2000