

Anti-CD74 antibody [SM4946] (STJA0004946)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H) 1:500-1:1000 IHC/IF (H) 1:500-1:1000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG1k
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION




Gene ID	972
Gene Symbol	CD74
UniProt ID	HG2A_HUMAN

ADDITIONAL INFORMATION




Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of CD74 at dilution of 1:1000



Immunohistochemistry analysis of paraffin-embedded human tonsil using CD74 at dilution of 1:1000



Immunofluorescent analysis of paraformaldehyde-fixed human tonsil using CD74 at dilution of 1:1000