

## Anti-Cytokeratin 5 antibody [SM4947] (STJA0004947)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES





Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H) 1:1000-1:2000 IHC/IF (H) 1:400-1:800
Formulation	PBS with 0.02% sodium azide, 100 µg/ml BSA and 50% glycerol.
Isotype	IgG2bk
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">110308</a>
Gene Symbol	<a href="#">Krt5</a>
UniProt ID	<a href="#">K2C5_MOUSE</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human esophagus using CK5 at dilution of 1:800

