

Anti-Recombinant-MYC-Tag antibody [SM4964] (STJA0004964)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/ICC/IF
Host / Source	Mouse
Reactivity	Myc-Tag

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:500-1:1000 ICC/IF (Species independent) 1:500-1:1000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of MYC Tag at dilution of 1:1000



Immunocytochemistry analysis of 4% paraformaldehyde fixed 293T cells using MYC Tag at dilution of 1:500