

## Anti-S100A6 antibody [SM5011] (STJA0005011)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC/IF             |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                                            |
|---------------------|------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                 |
| Conjugation         | Unconjugated                                                                                               |
| Purification        | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| Dilution Range      | IHC/IF (H) 1:200-1:500                                                                                     |
| Formulation         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                    |
| Storage Instruction | Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.                                    |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">6277</a>        |
| Gene Symbol | <a href="#">S100A6</a>      |
| UniProt ID  | <a href="#">S10A6_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

